


Benzodiazepines, Serum or Plasma, Quantitative

Order Name: **Benzodiaz S/P**
Test Number: 6980151
Revision Date: 01/04/2023

| TEST NAME | METHODOLOGY | LOINC CODE |
|--|--|------------|
| Benzodiazepines, Serum or Plasma, Quantitative | <u>Quantitative Liquid Chromatography/Tandem Mass Spectrometry</u> | |

| SPECIMEN REQUIREMENTS | | | | |
|-----------------------|--|---------------|----------------------------------|-----------------------|
| Specimen | Specimen Volume (mL) | Specimen Type | Specimen Container | Transport Environment |
| Preferred | 3 mL (1 mL) | Serum | Clot Activator (Red Top, No-Gel) | Refrigerated |
| Instructions | <p>Specimen Type: Red-top tube OR green-top (heparin) tube</p> <p>Specimen Storage: SUBMISSION/TRANSPORT (<3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen.</p> <p>Specimen Collection: Serum or plasma should be separated from cells within two hours of venipuncture. Send serum or plasma in a plastic transport tube.</p> <p>Special Instructions: Testing referred to MEDTOX Laboratories Inc TC19976</p> | | | |

| GENERAL INFORMATION | |
|---------------------|--|
| Expected TAT | 5-10 days |
| Notes | Labcorp Test Code: 791886 |
| CPT Code(s) | 80347 (Alt code: G0480) |
| Service Provided By |  labcorp Oklahoma, Inc. |