

**Lymphocyte Absolute CD4 Count**

Order Name: **T Helper**  
 Test Number: 3654600  
 Revision Date: 05/24/2023

TEST NAME	METHODOLOGY	LOINC CODE
Lymphocyte Absolute CD4 Count	Flow cytometry	

**SPECIMEN REQUIREMENTS**

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	<b>5 mL each (3 mL each)</b>	<b>Whole Blood</b>	<b>EDTA (Lavender) and ACD Sol A/B (Yellow Top)</b>	<b>Room Temperature</b>

Instructions	<p>Collect Whole blood in <b>Both ACD Solution A or B (Yellow Top)</b> and <b>EDTA (Lavender Top)</b>.            Keep whole blood, do not centrifuge or refrigerate.  <b>Deliver to the laboratory ASAP - Same Day as Collection!</b>  <b>Collect Monday through Wednesday only. Specimens must be received at Labcorp Oklahoma Inc. by Wednesday afternoon or collect on Monday.</b>  <b>Stability:</b> Room Temperature: Yellow-top (ACD-A) or (ACD-B) tube: Lavender-top (EDTA) tube: 1 day</p>
--------------	---

**GENERAL INFORMATION**

Testing Schedule	Testing schedule may vary.
Expected TAT	3-5 days
Clinical Use	Marker for level of immunocompetence. The T4 Lymphocytes reports out values for Lymphocyte Percent CD4 and Lymphocyte Absolute CD4
CPT Code(s)	86361
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.