


**Lymphocyte Absolute CD4 Count**

Order Name: **T Helper**  
Test Number: 3654600  
Revision Date: 05/24/2023

TEST NAME	METHODOLOGY	LOINC CODE
Lymphocyte Absolute CD4 Count	Flow cytometry	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	5 mL each (3 mL each)	Whole Blood	EDTA (Lavender) and ACD Sol A/B (Yellow Top)	Room Temperature
Instructions	Collect Whole blood in <b>Both ACD Solution A or B (Yellow Top)</b> and <b>EDTA (Lavender Top)</b> . Keep whole blood, do not centrifuge or refrigerate. <b>Deliver to the laboratory ASAP - Same Day as Collection!</b> <b>Collect Monday through Wednesday only. Specimens must be received at Labcorp Oklahoma Inc. by Wednesday afternoon or collect on Monday.</b> <b>Stability:</b> Room Temperature: Yellow-top (ACD-A) or (ACD-B) tube: Lavender-top (EDTA) tube: 1 day			

GENERAL INFORMATION	
Testing Schedule	Testing schedule may vary.
Expected TAT	3-5 days
Clinical Use	Marker for level of immunocompetence. The T4 Lymphocytes reports out values for Lymphocyte Percent CD4 and Lymphocyte Absolute CD4
CPT Code(s)	86361
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.