


**Ethanol Level Urine**

Order Name: **ETHANOL U**  
Test Number: 4314650  
Revision Date: 05/20/2021

TEST NAME	METHODOLOGY	LOINC CODE
Ethanol Urine	Enzymatic .....	42242-8

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	5 mL	Urine, Random	Sterile Urine container	Refrigerated
Instructions	Keep refrigerated. Specimen stability: Ambient 24 hours. Refrigerated 30 days.			

GENERAL INFORMATION	
Testing Schedule	Sun-Sat
Expected TAT	1-2 Days
Clinical Use	Detects presence of ethanol in the urine. For toxicity assessment, blood alcohol level is preferred. Urine ethanol may be positive due to the action of fungi or bacteria in urine without ingestion of ethanol, most commonly in diabetics.
CPT Code(s)	80307
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.