

## Pregabalin (Lyrica) Serum or Plasma

Order Name: Pregabalin S/P

Test Number: 6906325 Revision Date: 12/12/2022

TEST NAME			DDOLOGY	LOINC CODE
Pregabalin, Serum or Plasma		<u>lmmun</u>	passay (IA)	47414-8
SPECIMEN REQUIREM	MENTS			
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5 mL)	Serum	Clot Activator (Red Top, No-Gel)	Refrigerated
Alternate 1	1 mL (0.5 mL)	Plasma	Lithium Heparin (Dark Green Top / No-Gel)	Refrigerated
Instructions	Specimen Type: Red-top tube or green-top (heparin) tube. Gel-barrier tubes are NOT recommended.  Specimen Storage: SUBMISSION/TRANSPORT (<3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen.  Specimen Collection: Serum or plasma should be separated from cells within two hours of venipuncture. Submit serum or plasma in a plastic transport tube.			

GENERAL INFORMATION		
Expected TAT	2 - 7 days	
Notes	Labcorp Test Code: 808320	
CPT Code(s)	80366 (Alt code: G0480)	
Service Provided By	labcorp Oklahoma, Inc.	

Special Instructions: Testing referred to MEDTOX Laboratories Inc TC 237