

Leukotriene E4, 24-Hour Urine

Order Name: Leukot E4 Ur

Test Number: None
Revision Date: 10/25/2023

TEST NAME METHODOLOGY LOINC CODE

Leukotriene E4, 24-Hour Urine

Liquid Chromatography/Tandem Mass Spectrometry

SPECIMEN REQUIREMENTS					
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment	
Preferred	5 mL (2 mL)	Urine, 24-hour	24 hour Urine Container	Frozen	
Instructions	Order using New Test Misc: Use Labcorp test number 831775				
	Collection Instructions:				
	1. Start collection within a few hours of symptom onset and collect urine for 24 hours				
	2. No preservative preferred				
	3. Aliquot urine into a plastic v	ial and send frozen.			
	Stability: Room Temperature	24 hours; Refrigerated: 7 days;	Frozen: 28 days		

GENERAL INFORMATION		
Testing Schedule	Monday, Thursday	
Expected TAT	4-11 days after test set-up	
Clinical Use	Aiding in the evaluation of patients at-risk for mast cell activation syndrome (eg, systemic mastocytosis, IgE-mediated allergies, or aspirinexacerbated respiratory disease) using 24-hour urine collections	
Notes	Performing laboratory Mayo Medical Laboratories tc# TLTE4	
CPT Code(s)	82542	
Service Provided By	labcorp Oklahoma, Inc.	