

St. John Health System Lab Catalog

> Order Name: Leukot E4 Ur Test Number: None Revision Date: 10/25/2023

TEST NAME			METHODOLOGY	
Leukotriene E4, 24-Hour Urine			Liquid Chromatography/Tandem Mass Spectrometry	
SPECIMEN REQUIR	EMENTS			
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	5 mL (2 mL)	Urine, 24-hour	24 hour Urine Container	Frozen
Instructions	Uctions Order using New Test Misc: Use Labcorp test number 831775			
GENERAL INFORMA	Collection Instructions: Start collection within a few hours of symptom onset and collect urine for 24 hours No preservative preferred Aliquot urine into a plastic vial and send frozen. Stability: Room Temperature: 24 hours; Refrigerated: 7 days; Frozen: 28 days			
Testing Schedule	Monday, Thursday			
Expected TAT	4-11 days after test set-up			
Clinical Use	Aiding in the evaluation of patients at-risk for mast cell activation syndrome (eg, systemic mastocytosis, IgE-mediated allergies, or aspirin- exacerbated respiratory disease) using 24-hour urine collections			
Notes	Performing laboratory Mayo Medical Laboratories tc# TLTE4			

82542

labcorp Oklahoma, Inc.

CPT Code(s)

Service Provided By