


Porphyrins, Feces

Order Name: **Porph Feces**
Test Number: 5194403
Revision Date: 01/04/2025

| TEST NAME | METHODOLOGY | LOINC CODE |
|-------------------|--|------------|
| Porphyrins, Feces | High Performance Liquid Chromatography | |

| SPECIMEN REQUIREMENTS | | | | |
|-----------------------|---|------------------------|----------------------------|-----------------------|
| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment |
| Preferred | 100g | Timed Stool Collection | Sterile Screwtop Container | Frozen |
| Instructions | <p>Patient Preparation: Patient should refrain from ingesting red meat and aspirin-containing products (eg, Excedrin, Aspirin) for 3 days prior to specimen collection.</p> <p>Specimen: Entire collection (48, 72, or 96 hour) (min 10g) Stool in Sterile container Frozen, Note: 24-Hour collection is adequate if the collection volume is approximately 100 g. Length of collection period is required. Specimens smaller than 100 g may not provide interpretable results. Include a list of medications the patient is currently taking.</p> <p>Stability Requirements: Room temperature n/a, Refrigerated 14 days, Frozen 14 days</p> <p>Cause for Rejection: Specimen less than 10 grams.</p> <p>Special Instructions: Testing referred to Mayo Medical Lab TC#FQPPS</p> | | | |

| GENERAL INFORMATION | |
|---------------------|--|
| Expected TAT | 13 - 18 days |
| CPT Code(s) | 84126 |
| Service Provided By |  labcorp Oklahoma, Inc. |