

Endomysial Antibodies IgA, IgG

Order Name: **IgA IgG EMA**
Test Number: **5194823**
Revision Date: **12/12/2022**

TEST NAME	METHODOLOGY	LOINC CODE
Endomysial Antibodies IgA, IgG		

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5 mL)	Serum	Clot Activator (SST or Red No-Gel)	Room Temperature
Instructions	<p>Specimen Type: Gel-barrier tube or red-top tube</p> <p>Specimen Storage: Room temperature</p> <p>Specimen Collection: Serum should be separated from cells within two hours of venipuncture. Send serum in a plastic transport tube.</p> <p>Special Instructions: Testing referred to Beutner Labs BUFNY#505L</p> <p>Specimen Stability: Ambient: 5 days, Refrigerated : 5 days, Frozen: 1 year</p>			

GENERAL INFORMATION	
Expected TAT	7 days
Notes	Labcorp Test Code: 831723
CPT Code(s)	86256x2
Service Provided By	 labcorp Oklahoma, Inc.