


**Manganese, Urine**

Order Name: **Mang U Rand**  
Test Number: **5194547**  
Revision Date: **12/12/2022**

TEST NAME	METHODOLOGY	LOINC CODE
Manganese, Urine		

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Instructions	<p><b>Specimen Type:</b> Plastic urine container, no preservative</p> <p><b>Specimen Storage:</b> Maintain specimen at room temperature.</p> <p><b>Specimen Collection:</b> Not Available</p> <p><b>Special Instructions:</b> Do not use preservatives. Preservatives used for routine analysis may contain mercuric oxide (ie, Stabilur), which interferes with all metal testing. If both urinalysis and metal testing are ordered, please submit a separate urine specimen (containing no a</p> <p><b>Specimen Stability:</b> Ambient: 14 days, Refrigerated : 14 days, Frozen: 14 days</p>			

GENERAL INFORMATION	
Notes	Labcorp Test Code: 071597
CPT Code(s)	82570, 83785
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.