


**VMA, Random Urine**

Order Name: **VMA U Rand**  
Test Number: 5194550  
Revision Date: 12/12/2022

TEST NAME	METHODOLOGY	LOINC CODE
VMA, Random Urine		

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Instructions	<p><b>Specimen Type:</b> Plastic urine container</p> <p><b>Specimen Storage:</b> Vanillylmandelic acid is STABLE at room temperature unpreserved for seven days.</p> <p><b>Specimen Collection:</b> If specimen is to be stored for longer than seven days before analysis, use random urine transport tube.</p> <p><b>Special Instructions:</b> Use for patients where 24-hour urine collection is difficult (eg, pediatric patients). Twenty-four hour urine collection is preferable (see Vanillylmandelic Acid (VMA), 24-Hour Urine [004143]).</p> <p><b>Specimen Stability:</b> Ambient: 7 days, Refrigerated : 14 days, Frozen: 14 days</p>			

GENERAL INFORMATION	
Notes	Labcorp Test Code: 123208
CPT Code(s)	82570, 84585
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.