

Amphetamines, MS, Ur, MAT

Instructions

Order Name: Amph, MS, Ur, MAT

Test Number: 5195574 Revision Date: 01/09/2024

TEST NAME			METHODOLOGY	LOINC CODE
Amphetamines, MS, Ur, MAT			Mass Spectrometry (MS)	
SPECIMEN REQUIRE	MENTS			
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	30 mL (15 mL)	Urine, Random	Sterile Urine Container Preservative Free	Room Temperature

Specimen: 30mL(15mL) Random Urine Preservative free Urine Transport tube Room Temperature or Refrigerated. Cause for Rejection: Insufficient volume, no ID on container, urine from preservative tube.

SUBMISSION/TRANSPORT~(<3~days): Room~Temperature.~For~storage~beyond~three~days,~specimen~should~be~refrigerated~or~frozen.

Stability: Room Temperature 5 days, Refrigerated 14 days, Frozen 14 days.

GENERAL INFORMATION	
Expected TAT	5-8 days
CPT Code(s)	80359, 80324
Service Provided By	labcorp Oklahoma, Inc.