

Instructions

Order Name: Etoh Biomark Ur

Test Number: 5195577
Revision Date: 01/09/2024

**TEST NAME METHODOLOGY** LOINC CODE Ethanol Biomarkers, MS, Ur, MAT Mass Spectrometry (MS) **SPECIMEN REQUIREMENTS** Specimen Specimen Volume (min) Specimen Type Specimen Container Transport Environment 30 mL (15 mL) Sterile Urine Container Preferred Urine, Random **Room Temperature Preservative Free** 

Specimen: 30mL(15mL) Random Urine Preservative free Urine Transport tube Room Temperature or Refrigerated. Cause for Rejection: Insufficient volume, no ID on container, urine from preservative tube.

SUBMISSION/TRANSPORT (<3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen.

Stability: Room Temperature 5 days, Refrigerated 14 days, Frozen 14 days.

GENERAL INFORMATION	
Expected TAT	5-8 days
CPT Code(s)	80321
Service Provided By	Oklahoma, Inc.