


**ToxAssure Flex 22, Urine**

Order Name: **TA Flex 22 UR**  
Test Number: **5195611**  
Revision Date: **01/09/2024**

TEST NAME	METHODOLOGY	LOINC CODE
ToxAssure Flex 22, Urine	Immunoassay (IA) .....	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	30 mL (15 mL)	Urine, Random	Sterile Urine Container Preservative Free	Room Temperature
Instructions	Specimen: 30mL(15mL) Random Urine Preservative free Urine Transport tube Room Temperature or Refrigerated. Cause for Rejection: Insufficient volume, no ID on container, urine from preservative tube. SUBMISSION/TRANSPORT (<3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen. Stability: Room Temperature 14 days, Refrigerated 14 days, Frozen 6 weeks.			

GENERAL INFORMATION	
Expected TAT	5-8 days
CPT Code(s)	80347, 80357, 80354, 80307, 80362, 80348, 80368
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.