

St. John Health System Lab Catalog

> Order Name: Amph Confirm UR Test Number: 5195600 Revision Date: 01/11/2024

TEST NAME			METHODOLOGY	LOINC CODE
Amphetamines Confirmation, Ur			Liquid Chromatography/Tandem Mass Spectrometry	
SPECIMEN REQU	IREMENTS			
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	30 mL (15 mL)	Urine, Random	Sterile Urine Container Preservative Free	Room Temperature
Instructions	Specimen: 30mL(15mL) Random Urine Preservative free Urine Transport tube Room Temperature or Refrigerated. Cause for Rejection: Insufficient volume, no ID on container, urine from preservative tube. SUBMISSION/TRANSPORT (<3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen. Stability: Room Temperature 14 days, Refrigerated 14 days, Frozen 6 weeks.			

GENERAL INFORMATION	
Expected TAT	5-8 days
CPT Code(s)	80324, 80359
Service Provided By	Oklahoma, Inc.