


**Eosinophil Smear**

Order Name: **EOS SMEAR**  
Test Number: 0107800  
Revision Date: 11/12/2019

TEST NAME	METHODOLOGY	LOINC CODE
Eosinophil % For Eos Smear	Wright Stain	30378-4
Neutrophil % For Eos Smear	Wright Stain	13597-0
Epithelial For Eos Smear	Wright Stain	16822-9
Cellularity For Eos Smear	Microscopy	
Source For Smear For Eos	Prompt	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred		Nasal Smear	Glass Slides with Holder	Room Temperature
Alternate 1		Nasal swab	PNP Swab (Green Cap)	Room Temperature
Alternate 2		Nasal Wash	Sterile Screwtop Container	Room Temperature
Instructions	Nasal smear, Sputum or Bronchial Washings - Specimen stability of 24hrs - Keep Room Temperature.			

GENERAL INFORMATION	
Testing Schedule	Daily
Expected TAT	1 Day
CPT Code(s)	89190
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.