

Hemosiderin, Urine

Order Name: **HEMOSDRN U**
Test Number: 1000750
Revision Date: 11/09/2020

TEST NAME	METHODOLOGY	LOINC CODE
Hemosiderin, Urine	<u>Microscopy</u>	17783-2

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	10 mL (5 mL)	Urine, Random	Sterile Urine container	Refrigerated
Instructions	Early morning specimens preferred. Refrigerate or deliver to lab immediately. Collect in preservative-free container Monday through Friday only. Specimens must be received in lab section by Friday evening. Stability: Refrigerated 24hours.			

GENERAL INFORMATION	
Testing Schedule	Mon - Fri
Expected TAT	1 Day
CPT Code(s)	83070
Service Provided By	 labcorp Oklahoma, Inc.