

Chlamydia pneumoniae Antibody, IgG, IgM

 Order Name: **CHLAM AB**

Test Number: 5571650

Revision Date: 06/09/2003

| TEST NAME | METHODOLOGY | LOINC CODE |
|-----------------------------------|-------------------------------|------------|
| Chlamydia pneumoniae IgG Antibody | Indirect Fluorescent Antibody | 44981-9 |
| Chlamydia pneumoniae IgM Antibody | Indirect Fluorescent Antibody | 24004-4 |

SPECIMEN REQUIREMENTS

| | | | | |
|-----------|-----------------------|---------------|--------------------|-----------------------|
| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment |
| Preferred | 1 mL | Serum | Clot Activator SST | Refrigerated |

GENERAL INFORMATION

| | |
|---------------------|---|
| Testing Schedule | Mon - Fri |
| Expected TAT | 3 Days |
| Clinical Use | C. pneumoniae is an important respiratory pathogen as well as a risk factor for coronary heart disease. Approximately 40-50% of the adult population have positive IgG titers (>1:16) to C. pneumoniae. |
| CPT Code(s) | 86631/86632 |
| Service Provided By |  labcorp Oklahoma, Inc. |