

**Urea Nitrogen, Blood Post Dialysis**

 Order Name: **BUN POST**

Test Number: 2009750

Revision Date: 06/04/2021

TEST NAME	METHODOLOGY	LOINC CODE
Urea Nitrogen, Blood Post Dialysis	Urease/GLDH	11064-3

**SPECIMEN REQUIREMENTS**

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	<b>1 mL (0.5)</b>	<b>Plasma</b>	<b>Lithium Heparin PST (Light Green Top)</b>	<b>Refrigerated</b>
Alternate 1	<b>1 mL (0.5)</b>	<b>Serum</b>	<b>Clot Activator SST</b>	<b>Refrigerated</b>
Instructions	Stability: Ambient 7 days, Refrigerated 7 days, Frozen 1 year			

**GENERAL INFORMATION**

Testing Schedule	Daily
Expected TAT	1-2 days
Clinical Use	Post dialysis. Useful in the monitoring of dialysis.
CPT Code(s)	84520
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.