

Electrolytes, Urine Random

Order Name: **ELECT UR**
Test Number: 2012575
Revision Date: 05/20/2021

TEST NAME	METHODOLOGY	LOINC CODE
<u>Chloride Urine Random</u>	<u>Ion-Selective Electrode</u>	2078-4
<u>Potassium Urine Random</u>	<u>Ion-Selective Electrode</u>	2828-2
<u>Sodium Urine Random</u>	<u>Ion-Selective Electrode</u>	2955-3

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5)	Urine, Random	Sterile Urine container	Refrigerated
Instructions	Random urine collection. No preservative. Keep refrigerated. Specimen stability: Ambient 24 hours, Refrigerated 7 days, Frozen 7 days			

GENERAL INFORMATION	
Testing Schedule	Daily
Expected TAT	1-2 days
Clinical Use	Used to evaluate electrolyte balance, acid-base balance
CPT Code(s)	82436; 84300; 84133
Service Provided By	 labcorp Oklahoma, Inc.