


Urea Nitrogen, Blood (Arterial)

Order Name: **ART BUN**
Test Number: 2015650
Revision Date: 06/04/2021

| TEST NAME | METHODOLOGY | LOINC CODE |
|---------------------------------|----------------------|------------|
| Urea Nitrogen, Blood (Arterial) | Urease/GLDH | 12961-9 |

| SPECIMEN REQUIREMENTS | | | | |
|-----------------------|---|---------------|---------------------------------------|-----------------------|
| Specimen | Specimen Volume (mL) | Specimen Type | Specimen Container | Transport Environment |
| Preferred | 1 mL (0.5) | Plasma | Lithium Heparin PST (Light Green Top) | Refrigerated |
| Alternate 1 | 1 mL (0.5) | Serum | Clot Activator SST | Refrigerated |
| Instructions | Stability: Ambient 7 days, Refrigerated 7 days, Frozen 1 year | | | |

| GENERAL INFORMATION | |
|---------------------|--|
| Testing Schedule | Daily |
| Expected TAT | 1-2 days |
| Clinical Use | Useful in the diagnosis and treatment of certain renal diseases and metabolic disorders. |
| CPT Code(s) | 84520 |
| Service Provided By |  labcorp Oklahoma, Inc. |