


Creatinine, Serous Fluid

Order Name: **SRS CREAT**
Test Number: 2015850
Revision Date: 06/11/2003

| TEST NAME | METHODOLOGY | LOINC CODE |
|--------------------------|-----------------------------------------|------------|
| Creatinine, Serous Fluid | <u>Kinetic Alkaline Picrate (Jaffe)</u> | 12190-5 |

| SPECIMEN REQUIREMENTS | | | | |
|-----------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------|---------------|----------------------------|-----------------------|
| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment |
| Preferred | 1 mL (0.5) | Serous Fluid | Sterile Screwtop Container | Refrigerated |
| Instructions | Venous blood is often drawn simultaneously. Note fluid type on requisition and container. Specimen stability: Ambient 8 hours. Refrigerated 7 days. | | | |

| GENERAL INFORMATION | |
|---------------------|----------------------------------------------------------------------------------------------------------------------|
| Testing Schedule | Daily |
| Expected TAT | 1-2 days |
| Clinical Use | Useful for sensitive and quantitative measurement of renal function. |
| CPT Code(s) | 82570 |
| Service Provided By |  labcorp Oklahoma, Inc. |