


Creatinine, Urine Random

Order Name: **CREAT R U**
Test Number: 3000750
Revision Date: 11/04/2025

TEST NAME	METHODOLOGY	LOINC CODE
Creatinine, Urine Random	<u>Kinetic Alkaline Picrate (Jaffe)</u>	2161-8

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5)	Urine, Random	Sterile Urine container	Refrigerated
Instructions	Random urine collection. No preservative. Keep refrigerated. Specimen stability: Ambient 3 days, Refrigerated 3 days, Frozen 2 weeks			

GENERAL INFORMATION	
Testing Schedule	Daily
Expected TAT	1-2 days
Clinical Use	Useful in monitoring renal function, to follow possible progression of renal disease, adjust dosages of medications in which renal excretion is pivotal.
CPT Code(s)	82570
Service Provided By	 labcorp Oklahoma, Inc.