


Potassium Urine Random

Order Name: **POTAS R U**
Test Number: 3001350
Revision Date: 05/28/2021

TEST NAME	METHODOLOGY	LOINC CODE
Potassium Urine Random	<u>Ion-Selective Electrode</u>	2828-2

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5)	Urine, Random	Sterile Urine container	Refrigerated
Instructions	Random urine collection. No preservative. Keep refrigerated. Specimen stability: Ambient 24 hours, Refrigerated 2 months, Frozen 1 year			

GENERAL INFORMATION	
Testing Schedule	Daily
Expected TAT	1-2 days
Clinical Use	Useful in the management of kidney disease.
CPT Code(s)	84133
Service Provided By	 labcorp Oklahoma, Inc.