


Calcium Urine Random

Order Name: **CALC R U**
 Test Number: 3002250
 Revision Date: 04/30/2021

TEST NAME	METHODOLOGY	LOINC CODE
Calcium Urine Random	Arsenazo *****	17862-4

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5)	Urine, Random	Sterile Urine container	Refrigerated
Instructions	Random urine specimen. No preservative. Keep refrigerated. 2 - 25 ml of 6 N HCL or 3 - 10 ml of boric acid are acceptable preservatives if collecting with another test. Stability: Ambient 2 days, Refrigerated 4 days, Frozen 3 weeks.			

GENERAL INFORMATION	
Testing Schedule	Daily
Expected TAT	1-2 days
Clinical Use	Useful for identification of abnormal physiologic states causing excess or suppressed excretion of calcium, such as hyperparathyroidism, vitamin D abnormality, diseases that destroy bone, prostate cancer and drug treatment such as thiazide therapy.
CPT Code(s)	82310
Service Provided By	 labcorp Oklahoma, Inc.