

LDH Spinal FluidOrder Name: **CSF LDH**

Test Number: 3500700

Revision Date: 10/28/2011

| TEST NAME | METHODOLOGY | LOINC CODE |
|------------------|--------------------------|------------|
| LDH Spinal Fluid | Lactate - Pyruvate (NAD) | 60024-7 |

SPECIMEN REQUIREMENTS

| | | | | |
|--------------|--|----------------------------------|-----------------------------------|-----------------------|
| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment |
| Preferred | 0.5mL (0.1) | CSF (Cerebrospinal Fluid) | Sterile Screwtop Container | Refrigerated |
| Instructions | Patient should be informed, relaxed and properly positioned for lumbar puncture. Specimen stability: Ambient 6 hours. Refrigerated 24 hours. | | | |

GENERAL INFORMATION

| | |
|---------------------|---|
| Testing Schedule | Daily |
| Expected TAT | 1-2 days |
| CPT Code(s) | 83615 |
| Service Provided By |  labcorp Oklahoma, Inc. |