



## Acetaminophen Quantitative

Order Name: **ACETAMIN**  
Test Number: 4000050  
Revision Date: 04/26/2021

TEST NAME	METHODOLOGY	LOINC CODE
Acetaminophen Quantitative	Enzymatic	3298-7

### SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	<b>1 mL (0.5)</b>	<b>Plasma</b>	<b>Lithium Heparin PST (Light Green Top)</b>	<b>Refrigerated</b>
Alternate 1	<b>1 mL (0.5)</b>	<b>Serum</b>	<b>Clot Activator SST</b>	<b>Refrigerated</b>
Instructions	Recommended collection time is four hours after an oral dose. Specimen stability: Ambient 24 hrs, Refrigerated 14 days, Frozen 45 days			

### GENERAL INFORMATION

Testing Schedule	Daily
Expected TAT	1-2 days
Clinical Use	Useful for monitoring toxicity in overdose cases.
CPT Code(s)	80329
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.