

**Theophylline**

Order Name: **THEOPHYLLI**  
 Test Number: 4005050  
 Revision Date: 06/04/2021

TEST NAME	METHODOLOGY	LOINC CODE
Theophylline	Enzyme Immunoassay	4049-3

**SPECIMEN REQUIREMENTS**

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	<b>1 mL (0.5)</b>	<b>Plasma</b>	<b>Lithium Heparin PST (Light Green Top)</b>	<b>Refrigerated</b>
Alternate 1	<b>1 mL (0.5)</b>	<b>Serum</b>	<b>Clot Activator SST</b>	<b>Refrigerated</b>
Instructions	Peak: Draw 2-4 hours after oral dose. Stability: Ambient 8 hours, Refrigerated 3 months, Frozen 3 months.			

**GENERAL INFORMATION**

Testing Schedule	Daily
Expected TAT	1-2 days
Clinical Use	Useful for optimizing drug dosage and assessing toxicity.
CPT Code(s)	80198
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.