


**Tobramycin Peak**

Order Name: **TOBRA PEAK**  
Test Number: 4005600  
Revision Date: 09/20/2017

TEST NAME	METHODOLOGY	LOINC CODE
Tobramycin Peak	Enzyme Immunoassay	4057-6

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5)	Plasma	Lithium Heparin PST (Light Green Top)	Refrigerated
Alternate 1	1 mL (0.5)	Serum	Clot Activator SST	Refrigerated
Instructions	Peak: draw 30-60 minutes after 30 minutes infusion or 1 hour after IM dose. Stability: Ambient=8 hours; Refrigerated= 7 days			

GENERAL INFORMATION	
Testing Schedule	Daily
Expected TAT	1-2 days
Clinical Use	Useful for optimizing drug dosage and assessing toxicity.
CPT Code(s)	80200
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.