


**Benzodiazepine Screen Urine**

Order Name: **BENZO U**  
Test Number: 4301150  
Revision Date: 04/26/2021

TEST NAME	METHODOLOGY	LOINC CODE
Benzodiazepine Screen Urine	CEDIA	3390-2

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	10 mL (0.5)	Urine, Random	Sterile Urine container	Refrigerated
Instructions	Confirmations are done on positives upon request. Stability: Ambient 24 hours, Refrigerated 5 days, Frozen 1 month			

GENERAL INFORMATION	
Testing Schedule	Daily
Expected TAT	1-2 days
Clinical Use	Useful for monitoring toxicity in overdose cases.
CPT Code(s)	80307
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.