

Barbiturates Screen

Order Name: **BARB SC**

Test Number: 4301700

Revision Date: 04/26/2021

TEST NAME	METHODOLOGY	LOINC CODE
Barbiturates Screen	CEDIA	20421-4

SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5)	Serum	Clot Activator SST	Refrigerated
Alternate 1	1 mL (0.5)	Plasma	Lithium Heparin PST (Light Green Top)	Refrigerated
Instructions	Specimen stability: Ambient 24 hours, Refrigerated 5 days, Frozen 1 months			

GENERAL INFORMATION

Testing Schedule	Daily
Expected TAT	1-3 days
Clinical Use	Useful for monitoring toxicity in overdose cases.
CPT Code(s)	80307
Service Provided By	 labcorp Oklahoma, Inc.