

**Benzodiazepines Screen**

Order Name: **BENZ SC**  
Test Number: 4301800  
Revision Date: 04/26/2021

TEST NAME	METHODOLOGY	LOINC CODE
Benzodiazepines Screen	CEDIA	46976-7

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5)	Serum	Clot Activator SST	Refrigerated
Alternate 1	1 mL (0.5)	Plasma	Lithium Heparin PST (Light Green Top)	Refrigerated
Instructions	Specimen stability: Ambient 8 hours. Refrigerated 7 days, Frozen 1 month			

GENERAL INFORMATION	
Testing Schedule	Daily
Expected TAT	1-2 days
Clinical Use	Useful for monitoring toxicity in overdose cases.
CPT Code(s)	80307
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.