

Acetaminophen Screen Serum

Order Name: **ACETAM SC**
Test Number: 4302050
Revision Date: 04/26/2021

TEST NAME	METHODOLOGY	LOINC CODE
Acetaminophen Screen Serum	CEDIA	3297-9

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5)	Plasma	Lithium Heparin (Dark Green Top / No-Gel)	Refrigerated
Alternate 1	1 mL (0.5)	Serum	Clot Activator (Red Top, No-Gel)	Refrigerated
Instructions	Stability: Ambient 24 hrs, Refrigerated 14 days, Frozen 45 days			

GENERAL INFORMATION	
Testing Schedule	Daily
Expected TAT	1-2 days
Clinical Use	Useful for monitoring toxicity in overdose cases.
CPT Code(s)	80307
Service Provided By	 labcorp Oklahoma, Inc.