

**Flecainide (Tambocor)**

 Order Name: **TAMBOCOR**

Test Number: 3618200

Revision Date: 12/12/2022

TEST NAME	METHODOLOGY	LOINC CODE
<b>Flecainide (Tambocor)</b>	Liquid Chromatography/Tandem Mass Spectrometry	3638-4

**SPECIMEN REQUIREMENTS**

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	<b>1 mL (0.2 mL)</b>	<b>Serum</b>	<b>Clot Activator (Red Top, No-Gel)</b>	<b>Refrigerated</b>
Alternate 1	<b>1 mL (0.2 mL)</b>	<b>Plasma</b>	<b>Lithium Heparin (Dark Green Top / No-Gel)</b>	<b>Refrigerated</b>

Instructions	<b>Specimen Type:</b> Red-top (no additive) tube or green-top (heparin) tube; gel-barrier tubes not recommended <b>Specimen Storage:</b> Refrigerated <b>Specimen Collection:</b> Transfer separated serum or plasma to a plastic transport tube. Draw immediately prior to next dose unless otherwise instructed. <b>Specimen Stability:</b> Ambient: 3 days, Refrigerated : 14 days, Frozen: 14 days
--------------	---

**GENERAL INFORMATION**

Expected TAT	2 - 4 days working day.
Clinical Use	Monitoring the flecainide concentration is used to assure compliance and avoid toxicity of this cardiac drug used to treat ventricular tachcardia and premature cardiac contractions.
Notes	Labcorp Test Code: 085662
CPT Code(s)	80181
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.