


Beta-2, Transferrin BF

Order Name: **BETA-2 TRA**  
Test Number: 3656675  
Revision Date: 10/30/2025

TEST NAME	METHODOLOGY	LOINC CODE
Beta-2, Transferrin BF	Immunofixation Electrophoresis	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	2 mL	Body Fluid	Sterile Screwtop Container	Refrigerated
Instructions	<p>Aural or nasal fluid in sterile container; Do Not Freeze. Do Not send CSF specimens!</p> <p>Notes: 1 mL (Note: This volume Does NOT allow for repeat testing).</p> <p>Specimen Type: Sterile container</p> <p>Specimen Storage: Refrigerated</p> <p>Special Instructions: Testing referred to ARUP TC 50047</p> <p>Specimen Stability: Ambient: 2 weeks; Refrigerated: 2 weeks; Frozen: 2 weeks</p>			

GENERAL INFORMATION	
Expected TAT	4 - 8 days
Clinical Use	This test is to detect CSF in body fluids such as sinus or ear drainage.
CPT Code(s)	86335
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.