

**Lorazepam (Ativan)**

Order Name: **LORAZEPAM**  
Test Number: 3611975  
Revision Date: 12/20/2022

TEST NAME	METHODOLOGY	LOINC CODE
Lorazepam (Ativan)	<u>Quantitative Liquid Chromatography/Tandem Mass Spectrometry</u>	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	3 mL (0.6 mL)	Serum	Clot Activator SST	Room Temperature
Instructions	<p><b>Specimen Type:</b> Red-top tube OR green-top (heparin) tube. Gel-barrier tubes are NOT recommended.</p> <p><b>Specimen Storage:</b> SUBMISSION/TRANSPORT (&lt;3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen.</p> <p><b>Specimen Collection:</b> Serum or plasma should be separated from cells within two hours of venipuncture. Submit serum or plasma in a plastic transport tube.</p> <p><b>Special Instructions:</b> Testing referred to MedTox Laboratories Inc MX#64</p>			

GENERAL INFORMATION	
Expected TAT	7-10 Days
Notes	Labcorp Test Code: 808465
CPT Code(s)	80299
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.