


Chronic Urticaria PD-BATC

Order Name: **CU INDEX**  
Test Number: 5587555  
Revision Date: 06/03/2025

TEST NAME	METHODOLOGY	LOINC CODE
Chronic Urticaria PD-BATC	Flow cytometry .....	75023-2

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5 mL)	Serum	Clot Activator (SST or Red No-Gel)	Frozen
Instructions	<p><b>Specimen:</b> Serum 1 mL (0.5 mL Minimum Volume: 0.5 mL (Note: This volume does not allow for repeat testing.)</p> <p><b>Container:</b> Gel-barrier tube or red-top tube</p> <p><b>Collection:</b> <b>Allow to clot then Transfer separated serum to a plastic transport tube and Freeze.</b></p> <p><b>Storage Requirements:</b>Frozen</p> <p><b>Stability Requirements:</b> Room temperature 1 day, Refrigerated 14 days, Frozen 2 months (Freeze/thaw cycles Stable x3)</p> <p><b>Cause for Rejection:</b> Specimen that are inadequately labeled, lipemic, hemolyzed, cloudy or microbiologically contaminated.</p>			

GENERAL INFORMATION	
Expected TAT	7 - 10 days
Clinical Use	Chronic Urticaria - Pooled Donor Basophil Activation Testing (PD-BAT)
CPT Code(s)	86849
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.