

Cryofibrinogen

Order Name: **CRYOFIBRIN**
 Test Number: 5221675
 Revision Date: 12/12/2022

TEST NAME	METHODOLOGY	LOINC CODE
Cryofibrinogen	Precipitation	

SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	3mL (1mL)	Plasma	Sodium Citrate 3.2% (Blue Top)	See Instructions
Instructions	Specimen Type: Blue-top (sodium citrate) tube, lavender-top (EDTA) tube, or gray-top (sodium fluoride) tube Specimen Storage: Maintain specimen at room temperature. Specimen Collection: Specimen must be drawn in a prewarmed tube and kept at 37(degrees)C. Separate plasma from cells immediately and transfer the plasma into a clean transport tube. LABEL TRANSPORT TUBE "PREWARMED." Specimen Stability: Ambient: Not Available, Refrigerated : Not Available, Frozen: Not Available			

GENERAL INFORMATION

Expected TAT	3 - 5 days later.
Notes	Labcorp Test Code: 080317
CPT Code(s)	82585
Service Provided By	 labcorp Oklahoma, Inc.