


Ovarian Antibody Screen with reflex to Titer

Order Name: **OVARIAN AB**
Test Number: 5005750
Revision Date: 06/12/2009

TEST NAME	METHODOLOGY	LOINC CODE
Ovarian Antibody Screen with reflex to Titer	Immunofluorescence Assay (IFA)	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.3)	Serum	Clot Activator SST	Refrigerated
Instructions	Notes: 0.3 mL (Note: This volume Does not allow for repeat testing.) Specimen Type: Serum from a Gel Tube OR Serum from Red Top Tube Specimen Storage: Store frozen or refrigerated. Specimen Collection: Separate serum from cells within one hour. Transfer to a plastic transport tube before freezing. Special Instructions: Testing is performed at Esoterix Inc. ESECF#55 Specimen Stability: Ambient: 7 days, Refrigerated : 14 days, Frozen: 60 days			

GENERAL INFORMATION	
Expected TAT	7-14 days
Clinical Use	Ovarian Antibody is found in patients with premature ovarian failure, Addison's disease, and polyendocrinopathy syndromes.
CPT Code(s)	86256
Service Provided By	 labcorp Oklahoma, Inc.