

Helicobacter pylori CultureOrder Name: **C H PYLORI**

Test Number: 6002009

Revision Date: 09/21/2017

TEST NAME	METHODOLOGY	LOINC CODE
Helicobacter pylori Culture	Culture	

SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	Gastric biopsy		Sterile Screwtop Container	Refrigerated
Instructions	Place the tissue specimen and 1 mL of sterile saline in a sterile screw top container. Specimens in formalin or containers previously containing formalin are unacceptable.			

GENERAL INFORMATION

Testing Schedule	Daily
Expected TAT	10 Days
Clinical Use	Identifies Helicobacter pylori, causative agent of gastric ulcers
CPT Code(s)	87070
Service Provided By	 labcorp Oklahoma, Inc.