


Nasal Culture

Order Name: **C NASAL RT**
Test Number: 6002011
Revision Date: 09/27/2017

TEST NAME	METHODOLOGY	LOINC CODE
Nasal Culture	Culture	10353-1

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred		Nasal swab	Copan eSwab - White (Regular size)	Room Temperature
Alternate 1		Nasal swab	Copan eSwab - Green (Mini-tip)	Room Temperature
Alternate 2		Nasal swab	PNP Swab (Green Cap)	Room Temperature
Instructions	Swabs and aspirates should be transported at room temperature and setup within 2 days of collection for optimal recovery.			

GENERAL INFORMATION	
Testing Schedule	Daily
Expected TAT	3 Days
Clinical Use	Identifies upper respiratory pathogens
CPT Code(s)	87070
Service Provided By	 labcorp Oklahoma, Inc.