


**Venison IgE**

Order Name: **VENISON**  
Test Number: **5650145**  
Revision Date: **12/12/2022**

TEST NAME	METHODOLOGY	LOINC CODE
Venison IgE	Enzyme immunoassay (EIA)	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5 mL)	Serum	Clot Activator SST	Room Temperature
Instructions	<p><b>Specimen Type:</b> Serum transfer tube</p> <p><b>Specimen Storage:</b> Maintain specimen at room temperature</p> <p><b>Specimen Collection:</b> Serum must be separated from cells within 45 minutes of venipuncture. Send serum in a plastic transfer tube.</p> <p><b>Special Instructions:</b> Testing referred to Eurofins Viracor LLC TC 44710</p> <p><b>Specimen Stability:</b> Ambient: 4 weeks, Refrigerated : 4 weeks, Frozen: Not Available</p>			

GENERAL INFORMATION	
Expected TAT	6 - 10 days
Notes	Labcorp Test Code: 819091
CPT Code(s)	86003
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.