

## Venison IgE

Order Name: **VENISON**  
 Test Number: 5650145  
 Revision Date: 12/12/2022

TEST NAME	METHODOLOGY	LOINC CODE
Venison IgE	Enzyme immunoassay (EIA)	

### SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5 mL)	Serum	Clot Activator SST	Room Temperature
Instructions	<p><b>Specimen Type:</b> Serum transfer tube  <b>Specimen Storage:</b> Maintain specimen at room temperature  <b>Specimen Collection:</b> Serum must be separated from cells within 45 minutes of venipuncture. Send serum in a plastic transfer tube.  <b>Special Instructions:</b> Testing referred to Eurofins Viracor LLC TC 44710  <b>Specimen Stability:</b> Ambient: 4 weeks, Refrigerated : 4 weeks, Frozen: Not Available</p>			

### GENERAL INFORMATION

Expected TAT	6 - 10 days
Notes	Labcorp Test Code: 819091
CPT Code(s)	86003

Service Provided By

