

OK/TX General Allergy Panel

Order Name: **A OK/TX PN**
Test Number: **5606580**
Revision Date: **08/01/2023**

TEST NAME	METHODOLOGY	LOINC CODE
OK/TX General Allergy Panel	ImmunoCAP	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	10 mL (5 mL)	Serum	Clot Activator SST	Room Temperature
Instructions	Specimen Type: Red-top tube or gel-barrier tube, Separate Serum from Cells into a screwtop transport container. Stability Requirements: Room temperature 14 days, Refrigerated 14 days, Frozen 3 months. (Freeze/thaw cycles Stable x3)			

GENERAL INFORMATION

Expected TAT

3-5 Days

Notes

Alternaria alternata Allergen	Cottonwood Tree Allergen	Pecan, Hickory Tree Allergen
Aspergillus fumigatus Allergen	Dog Dander Allergen	Penicillin Chrysogenium Allergen
Bermuda Grass Allergen	Dust Mite (D. farinae) Allergen	Ragweed Common Short Allergen
Box Elder (Maple) Tree Allergen	Elm Tree Allergen	Rough Marshelder Allergen
Cat Dander Allergen	House Dust Mite Allergen	Sheep Sorrel Allergen
Cladosporium herbarum Allergen	Mountain Juniper Tree Allergen	Timothy Grass Allergen
Cockroach, German Allergen	Mulberry Tree White Allergen	White Ash Tree Allergen
Common Pigweed Allergen	Nettle Allergen	
Common Silver Birch Tree Allergen	Oak Tree Allergen	

Allergens performed at Labcorp Phoenix

CPT Code(s)

86003 x25

Service Provided By

labcorp
Oklahoma, Inc.