


**Follicular Lymphoma, IGH/BCL2, t(14;18) by FISH**

Order Name: **FOLLIC FSH**  
Test Number: 9116975  
Revision Date: 10/01/2022

TEST NAME	METHODOLOGY	LOINC CODE
Follicular Lymphoma, IGH/BCL2, t(14;18) by FISH	Fluorescence in Situ Hybridization .....	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	5 mL (1 mL)	Bone Marrow	Sodium Heparin (Green Top, No-Gel)	Room Temperature
Alternate 1	5 mL (3 mL)	Whole Blood	Sodium Heparin (Green Top, No-Gel)	Room Temperature
Instructions	Send specimen ASAP, Keep at room temperature! (DO NOT FREEZE). Frozen samples will be rejected.			

GENERAL INFORMATION	
Testing Schedule	Mon-Fri
Expected TAT	3-5 Days
Clinical Use	Useful to diagnose follicular lymphoma (> 75%) and some diffuse large B-cell lymphomas.
CPT Code(s)	88377
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.