

**Diuretics Screen, Urine**

Order Name: **DIURET S U**  
Test Number: 3632070  
Revision Date: 01/04/2023

TEST NAME	METHODOLOGY	LOINC CODE
Diuretics Screen, Urine	<a href="#">Liquid Chromatography/Tandem Mass Spectrometry</a>	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	10 mL (1.2 mL)	Urine, Random	Sterile Screwtop Container	Refrigerated
Instructions	<p><b>Specimen Type:</b> Plastic urine container without preservative</p> <p><b>Specimen Storage:</b> SUBMISSION/TRANSPORT (&lt;3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen.</p> <p><b>Special Instructions:</b> Testing referred to MedTox Laboratories Inc MX#413</p>			

GENERAL INFORMATION	
Expected TAT	10-14 Days
Notes	Labcorp Test Code: 808548
CPT Code(s)	80377
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.