


Acid Phosphate Total

Order Name: **ACID PHS T**
Test Number: 2000280
Revision Date: 12/12/2022

TEST NAME	METHODOLOGY	LOINC CODE
Acid Phosphate Total	Enzymatic	12173-1

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1.5 mL (0.8 mL)	Serum	Clot Activator (Red Top, No-Gel)	Frozen
Instructions	Allow specimen to clot completely. Transfer 1.5mL of serum into an aliquot tube and freeze. CRITICAL FROZEN. Separate specimens must be submitted when multiple tests are ordered. Unacceptable Conditions: Plasma. Non-frozen specimens. Hemolyzed specimens. Stability After separation from cells: Ambient: Unacceptable; Refrigerated: Unacceptable; Frozen: 1 month			

GENERAL INFORMATION	
Expected TAT	5 - 8 days
Clinical Use	This test is used in the evaluation of prostate cancer, Padgett's disease, hyperthyroidism, and Gaucher's disease.
Notes	Labcorp Test Code: 829642
CPT Code(s)	84060
Service Provided By	 labcorp Oklahoma, Inc.