

**Tomato Allergen IgG**

Order Name: **TOMATO IGG**  
Test Number: 3666200  
Revision Date: 12/14/2022

TEST NAME	METHODOLOGY	LOINC CODE
Tomato Allergen IgG		35547-9

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	0.3 mL (0.1 mL)	Serum	Clot Activator SST	Room Temperature
Instructions	<p><b>Specimen Type:</b> Red-top tube or gel-barrier tube</p> <p><b>Specimen Storage:</b> Room temperature</p> <p><b>Specimen Collection:</b> Not Available</p> <p><b>Special Instructions:</b> Values obtained with different assay methods should not be used interchangeably in serial testing. It is recommended that only one assay method be used consistently to monitor each patient's course of therapy.</p> <p><b>Specimen Stability:</b> Ambient: 14 days, Refrigerated : 14 days, Frozen: 3 months</p>			

GENERAL INFORMATION	
Expected TAT	2 - 4 days
Clinical Use	Tomato: <i>Lycopersicon lycopersicum</i>
Notes	Labcorp Test Code: 605051
CPT Code(s)	86001
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.