

Dandelion IgE

Order Name: **DANDELION**
Test Number: **5582300**
Revision Date: **12/09/2022**

TEST NAME	METHODOLOGY	LOINC CODE
Dandelion IgE	ImmunoCAP	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	0.2 mL (0.5 mL)	Serum	Clot Activator SST	Room Temperature
Instructions	Specimen Type: Red-top tube or gel-barrier tube Specimen Storage: Room temperature Specimen Collection: Not Available Specimen Stability: Ambient: Not Available, Refrigerated : Not Available, Frozen: Not Available			

GENERAL INFORMATION	
Expected TAT	2-4 Days
Clinical Use	Dandelion: <i>Taraxacum vulgare</i>
Notes	Labcorp Test Code: 602968
CPT Code(s)	86003
Service Provided By	 labcorp Oklahoma, Inc.