

21 Hydroxylase Antibody

Order Name: **21-Hydroxy Ab**
 Test Number: 3487055
 Revision Date: 11/14/2023

| TEST NAME | METHODOLOGY | LOINC CODE |
|-------------------------|-------------------|------------|
| 21 Hydroxylase Antibody | Qualitative ELISA | 85363-0 |

SPECIMEN REQUIREMENTS

| | | | | |
|--------------|--|---------------|---|-----------------------|
| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment |
| Preferred | 1.0mL (0.3mL) | Serum | Clot Activator SST | Frozen |
| Alternate 1 | 1.0mL (0.3mL) | Serum | Clot Activator (Red Top, No-Gel) | Frozen |
| Instructions | <p>Notes: 0.2 mL (Note: This volume Does NOT allow for repeat testing.)</p> <p>Specimen Type: Serum from Red Top Tube or Serum from a Gel Tube</p> <p>Specimen Storage: Frozen</p> <p>Specimen Collection: Separate serum from cells within one hour. Transfer to a plastic transport tube before freezing. To avoid delays in turnaround time when requesting multiple tests on frozen samples, please submit separate frozen specimens for each test requested.</p> <p>Special Instructions: Testing is performed at Esoterix Endocrinology Laboratory: ES#500036</p> <p>Specimen Stability: Ambient: 14 days, Refrigerated : 14 days, Frozen: 2 years</p> | | | |

GENERAL INFORMATION

| | |
|---------------------|--|
| Expected TAT | 8-11 days |
| Notes | Labcorp Test Code: 504805 |
| CPT Code(s) | 83516 |
| Service Provided By |  labcorp Oklahoma, Inc. |